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ROLF titanium frame is coated and does not need special care.

Avoid contact with products containing highly-concentrated solvents such as liquid cleaners, acetone, oils, vaseline or silicones. These products damage the frame's coating and the rubber ring in the FLEXLOCK hinge.

Never clean the frame with an ultrasonic cleaner.

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The titanium material and coating can show discoloration after some time. This is caused by ultraviolet radiation. Nevertheless, these changes in color are inevitable and form part of the individual patina. Furthermore, usage of cosmetics can cause color changes in the material and coating.

Please avoid exposing the frame to temperatures more than 40°C.

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There is the option to order any pair of glasses from the SKYLINE collection with customizable nose pieces and customizable temple tips.

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
Apart from two minor adjustments, the glazing is based on the tried and tested ROLF glazing system. There is no hole in the nose bridge. The glazing nylon thread is now inserted in each side of the lens. A small opening on the inside of each side of the the frame allows removal of the glazing thread and makes the lens exchange possible.

ROLF titanium frames must be glazed with high-index plastic lenses which have an edge thickness of a minimum of 2.2mm.


As the frame cannot be adjusted to the lenses, individually calculated lenses have to be used. When grinding the lenses, the correct size is obtained as follows: when inserting them into the frame a small tolerance is necessary. Subsequently it is entirely edged all around using a groove wheel with a width of 1,05mm to a groove depth of 0.5mm. If the lens is the correct size, it should remain in the frame without the filament; there should however be a small degree of tolerance. In case of doubt, it is recommended to test glazing with a trial lens.

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1. grind the lenses with the enclosed demo lenses.
2. insert lenses from the inside.
3. cut the glazing thread - length about 12cm - cut one side of the thread to a spike.

4. insert the spiky side into the opening first. Insert the glazing thread up to the nose bridge.
5. Push the thread out a little - cut it off - push it flush to the frame.

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The glazing thread can be carefully pushed out of the frame ' with the glazing tool which can be found in the provided First Aid Kit. Remove the glazing thread carefully, then remove the lenses.


**ATTENTION:** The tool is sharp and pointed. Watch your fingers! Also be careful that the coating of the spectacle frame is not damaged. In case of improper handling, ROLF does not assume warranty.

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Particular attention will be paid to our multidirectional FLEXLOCK hinges. A printed polyamide insert and a rubber ring connect the front of the glasses with the temples.

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In the flexible and screw less hinge joint, a highly elastic ring and a printed insert connect the front of the frame with the temples.

The FLEXLOCK hinge was awarded the Silmo d'Or 2017 in the category Engineering and Innovation.

*Offers an unprecedented  
level of flexibility!*

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
Optionally, the eyeglass frames can also be ordered with adjustable temple tips and adjustable nose pads.

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

Elegant-edel-schlicht-schick-dynamisch.

The models present themselves with a natural play of colors in raw, anthracite, gold, black and mocca. In addition, all versions are available in two finishes - Rough and Gentle.

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Titanium is a lightweight, extremely high quality and noble material and therefore ideal for the production of spectacles.

It provides stability in areas of highly sensitive technology, such as aviation and medicine.

Here, the material convinces with its resistance to external influences and its high skin tolerance. It owes its high durability and suitability for everyday use in the job as well as in the spare time to these characteristics.

Due to the 3D printing technology and the ROLF glazing system no solder joints are necessary and the front of the glasses does not have to be cut, which increases the stability of the spectacle frame and allows a reduced and timeless design.

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The Titanium Skyline Collection combines beautiful prescription glasses & sunglasses, and convinces with high-tech & futuristic design.

Through the harmonious transition from spectacle front to temple we create recognition to our previous collections.

Using design programs, our construction specialists and designers breathe life and elegance into the ROLF 3D-printed spectacles.

From the basic form to the technical innovations, such as hinges and glazing details, everything is designed and optimized in 3D programs.

After that, the whole thing has to pass our „design test“. Construction, fit and functionality are checked and approved for Titanium printing.

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Now the 3D printing starts. Powder and heat form the spectacle blanks.

The frames are then sanded, blasted, coated, tested and released for final processing by trained eyewear makers.

The 3D printing allows us to work almost independently of technical limits and designs. Compared to classic production, very complex geometries can be implemented.